

Analytical Balance (Digital)

Function:

- ✓ Electromagnetic force sensor
- ✓ LCD (White back light with black font)
- ✓ Tare function/Counting/Unit conversion
- ✓ Sensitivity/Speed set/Percentage
- ✓ Density function
- ✓ Accumulation function/peak hold/Upper and lower set



FA-HM Series

Model	FA2004HM	FA3004HM
Maximum Capacity:	230g	300g
Minimum Load:	10 mg	
Readability [d]:	0.1 mg	
Linearity:	±0.3g	
Stabilization time:	3S	
Display:	5.7 inches Color LCD (backlit), 4-Digit	
Key Pad	6 Keys	
Communication interface:	2×RS232,2×USB-A (interchangeable) USB-B, Wi-Fi@ (option)	
Power supply:	220 V AC 50/60 Hz	
Power Consumption:	Max. 10 W	
Operating temperature:	+10 to +40°C	
Weighing pan dimensions:	φ80/90 mm	
Packaging dimensions:	425 x 320 x 415 mm	
Weight:	7- 8 kg	
Standard Accessories:	Manual, AC adapter	