

**F-100MW-T Series Precision Electronic Weighing Scale****Main Technical Parameter:**

Mode	F-100MW-T
Range (g)	110g
Readability(g)	0.001g
Operating Temp	5°C-35°C
Pan Size	Ø90mm
Heat Source	Halogen Lamp
Temp.Sensor	PT 100
Temp. Range	40°C-199°C
Temp.Step	1°C
Moisture Range	0.00%-100.00%
Moisture Readability	0.01%
Dry Residual Range	100.00%-0.00%
Dry Residual Readability	0.01%
Temp Setting	40°C-199°C By 1°C Step
Time Setting	1-99min by 10 s
Number Of Storage	15
Number Of History	15
Heater	220V ±15% 50HZ 400W
Balance	Input: 220V 50HZ Output:9V
Host Dimension	360X490X350mm
The Chamber Height	22mm
Display	5 Inch screen
Calibration	External Calibration
Interface	RS232