

PP Fume Hood



LCD Display



PP Sink
resistant to strong
acid, alkali and
anti-corrosion



Waterproof Socket



Water&Gas
Remote Control

ADVANTAGE

1. It is safer to use anti-corrosive water tap.
2. Microprocessor control system, LCD display.
3. Made of porcelain white PP, resistant to strong acid, alkali and anti-corrosion.
4. Front window which is made of thick transparent toughened glass maximize light and visibility inside the fume hood, providing a bright and open working environment.

TECHNICAL PARAMETERS

Model		FSF1000P	FSF1200P	FSF1500P	FSF1800P
External Size (W*D*H)		1040*750*2200mm	1240*800*2200mm	1540*800*2200mm	1840*800*2200 mm
Internal Size (W*D*H)		820*520*872mm	1020*570*872mm	1320*570*872mm	1620*570*872 mm
Work Surface Height		700mm			
Max Opening		815mm			
Air Velocity		0.4~0.6m/s			
Noise		≤60dB			
Fluorescent Lamp		12W*1	30W*2	30W*2	30W*2
Blower		Built-in axial flow blower; Speed adjustable.			
Front Window		Manual, 5mm toughened glass, height adjustable.			
Power Supply		AC220V±10%, 50/60Hz; 110V±10%, 60Hz			
Consumption		330W	360W	360W	360W
Material	Main Body	Made of Porcelain white PP, thickness 8mm, resistant to strong acid, alkali and anti-corrosion			
	Work Table	Chemical resistant phenolic resin			
Standard Accessory		Fluorescent lamp, Water tap, Gas tap, Water Sink, Base cabinet			
		Waterproof socket*2 4 meters PVC duct, Diameter: 250mm			
Optional Accessory		PP work table, epoxy resin board or ceramic board. Outer centrifugal blower, resistant to strong acid, alkali.			
Gross Weight		160kg	198kg	225kg	259kg
Package Size (W*D*H)	Main Body	1110*900*1870mm	1310*885*1575mm	1610*885*1575mm	1910*885*1575mm
	Base Cabinet	1110*790*950mm	1310*840*950mm	1610*840*950mm	1910*840*1575mm