

## Digital Medical Infrared Thermometers

### Characteristics:

"body or objects surface temperature testing

one button operation

32\*2 memory groups for dual modes

3 colors back light

°C and °F switchable"



### Main Technical Parameters:

<b>Model</b>	JZK-601
<b>Memory groups</b>	32*2 groups
<b>Product Dimensions</b>	14.3* 9.3* 3.8 CM
<b>Business type</b>	Manufacturer
<b>Storage temperature</b>	-20-50°C
<b>Storage Humidity</b>	15%-90%RH
<b>Operating temperature</b>	5-40°C
<b>Operating Humidity</b>	40%-85%RH
<b>Display</b>	LCD crystal digital display
<b>Measuring method</b>	Infrared radiation temperature measurement
<b>Measurement range</b>	32~42.5°C fro human body, 0-100°C for objects
<b>Resolution</b>	0.1 degrees Celsius (0.1 degrees Fahrenheit)
<b>Accuracy</b>	0-39.9kpa(0-299 mmhg)
<b>Voltage</b>	3 V
<b>Power consumption</b>	< 3 mW
<b>Power supply</b>	2* AA battery
<b>Selling Units</b>	Single item
<b>Package Size</b>	14X20X6.5 cm
<b>Weight</b>	0.3kg

