

FAITHFUL[®] FAITHFUL INSTRUMENT (Hebei) CO.,Ltd

- Melting point is the temperature of something turning to liquid from solid. Testing it is the main method to detect some characters such as purity etc. It is suitable for testing Melting points of drug, spice and dye etc.

Applicable standards:

- Corporate standard(Thickness tester)

Main Features:

- Standard mercury thermometers measure temperature.
- Manual knob to change heater voltage
- Pointer voltmeter shows heater voltage
- Observe the melting of samples through a magnifying glass
- The melting point of three samples can be tested simultaneously.

Specifications:

- Range of melting point (50~300) °C
- Accuracy ±1.0°C
- Usage temperature (15~35) °C
- Rate (1.0、 2.0、 4.0、 6.0) °C / min
- Accuracy of Rate ≤10%
- Power of the heater 100W
- Power 220V / 50Hz /100W
- Dimension (190*120*200) mm³

