

SHZ-D III Circulating Water Vacuum Pump

1. Equipped with circulating device, vacuum degree is stable, appearance is beauty, small and exquisite.
2. This product has double tapping, can be single used or parallel used with two vacuum meters.
3. Main engines is made from stainless steel material, box is made form engineering plastics. Anti-corrosion, no pollution, low noise and convenient to move.



SHZ-D III (A)
Anti-Corrosion



SHZ-D III (P)
Teflon



SHZ-D III (S)
Stainless Steel

Model	Power	Voltage	Flow	Lift	Machine Material	Tapping	Single Tap Sucking Rate	Max Vacuum Degree	Water Storage Tank Volume	Dimension (L*W*H)	Weight
SHZ-D III (A)	180W	220V	60L/MIN	8M	Anti-Corrosion	2 pcs	10L/MIN	0.098MPA	15L	400*290*440	15KG
SHZ-D III (P)	180W	220V	60L/MIN	8M	Teflon	2 pcs	10L/MIN		15L	400*290*440	15KG
SHZ-D III (S)	180W	220V	60L/MIN	8M	Stainless Steel	2 pcs	10L/MIN		15L	400*290*440	15KG