

## Top Loading Balance (Digital)

### Function:

- ✓ Fashionable appearance
- ✓ Large size stainless steel pan
- ✓ LCD (White back light with black font)
- ✓ AC and DC power supply
- ✓ Tare function/Counting/Unit conversion (g/ct/oz)
- ✓ Min weighing set
- ✓ Overload alarm/Level indicator
- ✓ Option: Dual display/Interface/Printer/Dust cover



<b>FA-TLB Series</b>		
Model	FA1200TLB	FA2100TLB
Maximum Capacity:	1200g	2100g
Minimum unit weight	0.1g	
Repeatability:	0.001g	
Linearity:	±0.003g	
Sensitivity Drift	±20ppm/°C (10°C~30°C /50°F~86°F)	
Display:	7 segment LCD with backlight (character height 16mm)	
Display Refresh	10 times/second	
Load	0.02g	
Communication interface:	2×RS232,2×USB-A (interchangeable) USB-B, Wi-Fi@ (option)	
Power supply:	220 V AC 50/60 Hz	
Power Consumption:	Max.10 W	
Operating temperature:	+10 to +40°C	
Pan size (W) × (D):	133X170 mm	160 × 160 mm
Weight:	7- 8 kg	
Standard Accessories:	Manual, AC adapter	