

### Halogen Moisture Analyzer



#### Characteristic:

- Halogen light heating
- Stainless steel chamber
- Temperature and time can be set
- The percentage of moisture content
- The percentage of dry residual
- LCD with backlight
- Store historical setting

#### Main Technical Parameter:

Model	F-110MW	F-105MW	F-102MW	F-100MW	F-100MW-T
Capacity(g)	110g/10mg	110g/5mg	110g/2mg	110g/1mg	110g/1mg
Sensor	Load Cell			Load Cell	
Pan size	Ø90mm				
Read ability	10mg	5mg	2mg	1mg	1mg
Operation temp. range	5-35°C				
Heat source	HALOGEN LAMP-1			HALOGEN LAMP-2	
Temp sensor	PT-100				
Temp range	40-199°C				
Temp step	1°C				
Moisture range	0.00%-100.00%				
Moisture readability	0.2%	0.10%	0.04%	0.01%	0.01%
Dry range	100.00%-0.00%			100.00%-0.00%	100.00%-0.00%
Dry readability	0.2%	0.10%	0.04%	0.01%	0.01%
ATRO	100%...999%				100%...999%
ATRO	0%...999%				0%...999%
Number of storage	15				20
Packing size	490*360*350mm				
G. weight	7.5kg				