

HGS-5000 Series Glass Bead Sterilizer

Features:

Glass Bead Sterilizers are a Quick, Easy and Accurate alternative to traditional methods of sterilization, and compact enough to be placed on any Laminar Air Flow Work benches or any other workside tables in a clean room atmosphere.



HGS-5000A



HGS-5000B

Specification

Model	HGS-5000A	HGS-5000B
Crucible size	Φ40mm x 80mm	Φ40mm x 140mm
Weight of Glass Bead	150g	300g
Temperature Range	100°C~300°C	100°C~300°C
Temp. Control Accuracy	≤5°C	≤5°C
Heating Time (40°C~300°C)	≤25min	≤25min
Power	120W	250W
Weight	2kg	2.5kg
Dimensions (mm)	130 x 145 x 155 mm	130 x 145 x 215 mm

Order Information

No.	Code	Description
01	AS-07070-00	GBS-5000A, AC200-AC240V/AC100-AC120V, 50/60Hz
02	AS-07080-00	GBS-5000B, AC200-AC240V/AC100-AC120V, 50/60Hz
03	AS-07071-01	150g Glass beads (Φ3mm)
04	AS-07071-02	300g Glass beads (Φ3mm)